



## 24-6106

Kölbel Blade, 36 x 85mm, each



- Original 36 mm width
- Interchangeable Kölbel blades fit on any Shoulder Frame
- Good for larger patients
- Components of Standard and Deluxe Shoulder sets (24-6001 + 24-8001)
- 85 mm deep



Kölbel blades are distinguished by their peristronic countenance that provides a gentle uplift to the soft tissue and muscle it spreads. Ballooning blades extend a full 20 mm from the central axis of their hub before beginning the return sweep to their axis of origin. Excellent visualization and secure retention of the retractor in the wound are assured while shoulder musculature can often be more effectively retracted by pairing blades of differing lengths.