



24-6110

Kolbel Retractor Set w/6 blade



- Features popular and versatile Kölbel frame
- Includes pair of 36x36 mm deep blades
- Includes pair of 36x53 mm deep blades
- Includes pair of 36x68 mm deep blades
- All items are components of Deluxe and Optimal Shoulder Sets

One-handed placement, positioning and operation make the Kölbel Retractor Sets well suited for a broad range of open and mini-open shoulder repair procedures. Spring loaded ball bearings set into each arm of the retractor frame hold blades while allowing limitless free rotation.