



24-6368

MiniOpen Set, w/3 + 6 blades



- Modified Mini-Open retractor frame features a wingnut as a center screw
- Attach one of three included middle arms to the frame using the wingnut
- Middle blades depths are 36 mm, 53 mm and 68 mm.
- Includes one pair of 36x36 mm blades
- Includes one pair of 36x53 mm blades
- Includes one pair of 36x68 mm blades
- All items are components of Optimal Shoulder set (24-8500)

Modified Mini-Open Retractor Set
(frame, 3 middle blades & 6 arm
blades, 36 - 68 mm)