

# Buxton BioMedical




Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 345-211-50

### Bo-Peep Hook, 5.0 mm



- 220 mm overall length
- 130 mm shaft
- Slim, knurled handle
- 5 mm wide blunt hook

Lil' Bo-Peep Suture Retrieving Hook

[http://buxtonbio.com/bopeep\\_hook\\_5\\_0\\_mm\\_640](http://buxtonbio.com/bopeep_hook_5_0_mm_640)