

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

52-5255

Sq Hndl Hook,shrp,3-prong,7mm



- Overall length: 175 mm
- Shaft length: 75 mm
- 7 mm reach
- Triple Prong
- Sharp

Square Handle Hook, sharp, 3-prong,
7mm

http://buxtonbio.com/sq_hndl_hookshrp3prong7mm_612