

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

77-4205

Dahmen-Lexer Gouge, 5mm wide



- 225 mm overall length
- 112 mm shaft
- Autoclavable, linen-phenolic handle
- 5 mm wide gouge

Dahmen-Lexer Gouge, 5mm wide

http://buxtonbio.com/dahmenlexer_gouge_5mm_wide_556