

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

83-8404

W'l-Blakesley Rngur, dn 30°, 4mm



- Standard IVD
- Tear Drop Cup with Fenestration
- Down 30° Cup Angle
- 200 mm Shaft Length

Weil-Blakesley Rongeur, down 30°,
4.0mm

http://buxtonbio.com/wlblakesley_rngurdn_304mm_859