

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 833-161-15

IVD Rngr, L-hndl, 160x1.5, str



- Bayonet Shaft
- Offset Left
- Straight Cup Angle
- Shaft Length 160

IVD Rongeur, left handled, 160mm  
shaft, 1.5mm cup straight

[http://buxtonbio.com/ivd\\_rngr\\_lhndl\\_160x1\\_5\\_str\\_929](http://buxtonbio.com/ivd_rngr_lhndl_160x1_5_str_929)