

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 833-164-40

IVD Rngr, L-hndl, 160x4mm,dn30



- Specialty IVD
- Specialty Bayonet
- Offset Left
- Down 30° Cup Angle
- Shaft Length 160

IVD Rongeur, left handled, 160mm  
shaft, 4mm cup, down 30°

[http://buxtonbio.com/ivd\\_rngr\\_lhndl\\_160x4mmdn30\\_935](http://buxtonbio.com/ivd_rngr_lhndl_160x4mmdn30_935)