

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-6506

Caspar IVD Rngr, Cer, str, 140x6mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Straight Cup Angle
- 140 mm Shaft Length

Caspar IVD Rongeur, Ceramic coated,  
straight, 140x6mm, fenestrated cups

[http://buxtonbio.com/casparivd\\_rngrcerstr140x6mm\\_899](http://buxtonbio.com/casparivd_rngrcerstr140x6mm_899)