

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-6515

Caspar IVD Rngr, Cer, up 30° 140x5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Up 30° Cup Angle
- 140 mm Shaft Length

Caspar IVD Rongeur, Ceramic coated,  
up 30°, 140x5mm, fenestrated cups

[http://buxtonbio.com/casparivdrngrcerup30140x5mm\\_902](http://buxtonbio.com/casparivdrngrcerup30140x5mm_902)