

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-7481

Std IVD Rngr, Cer, dn 30° 300x3mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Down 30° Cup Angle
- 300 mm Shaft Length

Standard IVD Rongeur, Ceramic coated, down 30°, 300 x 3 mm

http://buxtonbio.com/std_ivd_rngrcerdn_30300x3mm_820