

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-7551

Weil-Bsly Rngr, Cer, str, 300x3mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Straight Cup Angle
- 300 mm Shaft Length

Weil-Blakesley Rongeurs, Ceramic coated, straight, 300 x 3 mm

[http://buxtonbio.com/weilbsly\\_rngrcerstr300x3mm\\_887](http://buxtonbio.com/weilbsly_rngrcerstr300x3mm_887)