

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-7563

Weil-Bsly Rngr, Cer, up30° 300x5



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Up 30° Cup Angle
- 300 mm Shaft Length

Weil-Blakesley Rongeurs, Ceramic coated, up 30°, 300 x 5 mm

[http://buxtonbio.com/weilbsly\\_rngrcerup30\\_300x5\\_889](http://buxtonbio.com/weilbsly_rngrcerup30_300x5_889)