

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-7573

Weil-Bsly Rngr, Cer, up45° 300x5



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Tear Drop Cup with Fenestration
- Up 45° Cup Angle
- 300 mm Shaft Length

Weil-Blakesley Rongeurs, Ceramic coated, up 45°, 300 x 5 mm

[http://buxtonbio.com/weilbsly\\_rngrcerup45\\_300x5\\_891](http://buxtonbio.com/weilbsly_rngrcerup45_300x5_891)