

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-7683

Caspar IVD Rngr, Cer, dn, 300x3mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 300 mm Shaft Length

Caspar IVD Rongeur, Ceramic coated,  
down angled, 300x3mm, fenestrated  
cups

[http://buxtonbio.com/caspar\\_ivd\\_rngrcerdn300x3mm\\_10109](http://buxtonbio.com/caspar_ivd_rngrcerdn300x3mm_10109)