

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

837-7684

Caspar IVD Rngr, Cer, dn, 300x4mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Double Fenestrated Cup with Teeth
- Down 30° Cup Angle
- 300 mm Shaft Length

Caspar IVD Rongeur, Ceramic coated,
down angled, 300x4mm, fenestrated
cups

http://buxtonbio.com/caspar_ivd_rngrcerdn300x4mm_10109