

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 837-7892

Peapod Rngr, Cer, down, 180x1.5mm



- Standard IVD
- Ceramic Black Finish
- Finger Ring Handle
- Shorter Cup
- Down 30° Cup Angle
- 180 mm Shaft Length

<div  
class="newsArtCenterBold">Peapod  
IVD Cup Rongeur</div>

[http://buxtonbio.com/peapod\\_rngrcerdown180x1\\_5mm\\_796](http://buxtonbio.com/peapod_rngrcerdown180x1_5mm_796)