

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 838-132-15

IVD Rngr, Cer, R-Hd 130x1.5, up30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Right
- Up 30° Cup Angle
- Shaft Length 130

IVD Rongeur, Ceramic coated, right handled, 130mm shaft, 1.5mm cup, up 30°

[http://buxtonbio.com/ivd\\_rngrcerrhd\\_130x1\\_5up30\\_970](http://buxtonbio.com/ivd_rngrcerrhd_130x1_5up30_970)