

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 838-134-40

IVD Rngr,Cer,R-Hndl 130x4,dn30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Right
- Down 30° Cup Angle
- Shaft Length 130

IVD Rongeur, Ceramic coated, right handled, 130mm shaft, 4mm cup, down 30°

[http://buxtonbio.com/ivd\\_rngrcerrhdl\\_130x4dn30\\_973](http://buxtonbio.com/ivd_rngrcerrhdl_130x4dn30_973)