

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 838-181-40

IVD Rngr, Cer, R-hnd, 180x4mm, str



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Right
- Straight Cup Angle
- Shaft Length 180

IVD Rongeur, Ceramic coated, right handled, 180mm shaft, 4mm cup straight

[http://buxtonbio.com/ivd\\_rngrcerrhnd180x4mmstr\\_10093](http://buxtonbio.com/ivd_rngrcerrhnd180x4mmstr_10093)