

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

839-112-15

IVD Rngr Cer L-hd 110x1.5,up30



- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Left
- Up 30° Cup Angle
- Shaft Length 110

IVD Rongeur, Ceramic coated, left handled, 110mm shaft, 1.5mm cup, up 30°

http://buxtonbio.com/ivd_rngr_cer_lhd_110x1_5up30_959