

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 839-134-20

IVD Rngr Cer L-hndl 130x2,dn30



Photo  
Pending

- Specialty IVD
- Ceramic Black Finish
- Specialty Bayonet
- Offset Left
- Down 30° Cup Angle
- Shaft Length 130

IVD Rongeur, Ceramic coated, left handled, 130mm shaft, 2mm cup, down 30°

[http://buxtonbio.com/ivd\\_rngr\\_cer\\_lhndl\\_130x2dn30\\_963](http://buxtonbio.com/ivd_rngr_cer_lhndl_130x2dn30_963)