

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

84-4031

Frs-Sm Punch, 180x3mm, 90°up



- Original Construction
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling None
- Ferris-Smith Handle
- 180 mm Shaft Length

Ferris-Smith Punch, 180x3mm, 90° up-angled

http://buxtonbio.com/frssm_punch_180x3mm90up_1036