

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

84-4244

Frs-Sm Punch, 200x4mm, 40°dwn



- Original Construction
- Down 40° Angled Foot Plate
- Lumbar; Regular Bite
- Ferris-Smith Handle
- 200 mm Shaft Length

Ferris-Smith Punch, 200x4mm, 40°
down-angled

http://buxtonbio.com/frssm_punch_200x4mm40dwn_10084