

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 84-5414

Punch, 230x1mm, 40° dwn, TFP



- Original Construction
- Down 40° Angled Foot Plate
- Cervical; Regular Bite
- 230 mm Shaft Length

Sm-Hndl Punch, 230x1mm, 40° down-angled, TFP

[http://buxtonbio.com/punch\\_230x1mm40dwntfp\\_1089](http://buxtonbio.com/punch_230x1mm40dwntfp_1089)