

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

84-6462

Sm-Hndl Punch, 230x6mm, 40° up



- Original Construction
- Up 40° Angled Foot Plate
- Lumbar; Deep Bite
- Chip Handling Chamber
- 230 mm Shaft Length

Chmbr Sm-Handle Punch, 230x6mm,
40° up-angled

http://buxtonbio.com/smhndl_punch_230x6mm40up_1025