

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-2461

LrgHnd Pnch, 230x6mm, 90° up, Ejec



- Larger 16 mm opening
- Ejector tip clears bone chips with each bite
- Up 90° Angled Foot Plate
- 230 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, 230x6mm, 90° up-angled, Ejec

[http://buxtonbio.com/lrg\\_hnd\\_pnch230x6mm90upejec\\_1059](http://buxtonbio.com/lrg_hnd_pnch230x6mm90upejec_1059)