

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-2522

LrgHnd Pnch, 250x2mm, 40° up, Ejec



- Larger 16 mm opening
- Ejector tip clears bone chips with each bite
- Up 40° Angled Foot Plate
- 250 mm Shaft Length
- 2 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, 250x2mm, 40° up-angled, Ejec

http://buxtonbio.com/lrg_hnd_pnch250x2mm40upejec_1016