

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 843-7542

Bay-Kers'nUp, RR, 40°, 110x4mm



- Satin
- Original
- Offset Up
- Shaft Length 110 mm
- Chip Clearing Ejector
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, 110 mm shaft, offset up, Ramrod ejector, 40° x4mm bite

[http://buxtonbio.com/baykersnuprr40110x4mm\\_1475](http://buxtonbio.com/baykersnuprr40110x4mm_1475)