

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-7722

Bay-Kers'nUp, RR, 40°, 180x2mm



- Satin
- Original
- Offset Up
- Shaft Length 180mm
- Chip Clearing Ejector
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, 180 mm shaft, offset up, Ramrod ejector, 40° x2mm bite

http://buxtonbio.com/baykersnuprr40180x2mm_1491