

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

843-8732

Bay-Ker, L-hnd, RR, 180x3, 40°



Photo
Pending

- Satin
- Original
- Offset Left Handed
- Shaft Length 180mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bayonet Kerrison rongeur, 180 mm shaft, left-handed, Ramrod ejector, 40° x3mm bite

http://buxtonbio.com/baykerlhndrr180x340_1519