

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

845-5164

Foramin. Rngur, lite cv, 4mm, tfp



- Foraminotomy punch with light "shallow" curve
- 200 mm shaft
- Satin Finish
- Original Construction
- Thin Footplate
- 4 mm bite

Foraminotomy Rongeurs, light curve,
4mm, thin footplate

http://buxtonbio.com/foramin_rngur_lite_cv4mmtfp_1449