

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2022

LrgHnd Punch, Cer, 180x2mm, 40°up



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 40° Angled Foot Plate
- 180 mm Shaft Length
- 2 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-HndI Punch, Cer, 180x2mm, 40°  
up-angle

[http://buxtonbio.com/lrghnd\\_punchcer180x2mm40up\\_1218](http://buxtonbio.com/lrghnd_punchcer180x2mm40up_1218)