

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2032

LrgHnd Punch, Cer, 180x3mm, 40°up

Photo
Pending

- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 40° Angled Foot Plate
- 180 mm Shaft Length
- 3 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 180x3mm, 40°
up-angle

http://buxtonbio.com/lrghnd_punchcer180x3mm40up_1218