

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2033

LrgHnd Pnch, Cer, 180x3mm, 90°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 90° Angled Foot Plate
- 180 mm Shaft Length
- 3 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 180x3mm, 90°
down-angle

http://buxtonbio.com/lrgHnd_pnchcer180x3mm90dwn_1330