

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2253

LrgHnd Pnch, Cer, 200x5mm, 90°dwn

Photo
Pending

- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 90° Angled Foot Plate
- 200 mm Shaft Length
- 5 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 200x5mm, 90°
down-angle

http://buxtonbio.com/lrghnd_pnchcer200x5mm90dwn_1331