

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2464

LrgHnd Pnch, Cer, 230x6mm, 40°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 40° Angled Foot Plate
- 230 mm Shaft Length
- 6 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 230x6mm, 40°  
down-angle

[http://buxtonbio.com/lrgwnd\\_pnchcer230x6mm40dwn\\_1310](http://buxtonbio.com/lrgwnd_pnchcer230x6mm40dwn_1310)