

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2532

Sm-Hdl Punch, Cer, 250x3mm, 40° up



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 40° Angled Foot Plate
- 250 mm Shaft Length
- 3 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Sm-Hndl Punch, Ceramic, 250x3mm,  
40° up-angle

[http://buxtonbio.com/smhdl\\_punchcer250x3mm40up\\_1221](http://buxtonbio.com/smhdl_punchcer250x3mm40up_1221)