

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2542

Sm-Hdl Punch, Cer, 250x4mm, 40° up



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 40° Angled Foot Plate
- 250 mm Shaft Length
- 4 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Sm-Hndl Punch, Ceramic, 250x4mm,
40° up-angle

http://buxtonbio.com/smhdl_punchcer250x4mm40up_1221