

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 847-2544

Punch, Cer, 250x4mm, 40°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 40° Angled Foot Plate
- 250 mm Shaft Length
- 4 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Sm-Hndl Punch, Cer, 250x4mm, 40°  
down-angle

[http://buxtonbio.com/punchcer250x4mm40dwn\\_1311](http://buxtonbio.com/punchcer250x4mm40dwn_1311)