

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2553

Punch, Cer, 250x5mm, 90°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 90° Angled Foot Plate
- 250 mm Shaft Length
- 5 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Sm-Hndl Punch, Cer, 250x5mm, 90°
down-angle

http://buxtonbio.com/punchcer250x5mm90dwn_1333