

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2822

Cer Ker Rngur,lg,40°up,280x2mm



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Up 40° Angled Foot Plate
- 180 mm Shaft Length
- 2 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands
- 180 mm Shaft Length

Cer Ker Rngur,lg,40°up,280x2mm

http://buxtonbio.com/cer_ker_rngurlg40up280x2mm_1218