

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-2944

LrgHnd Pnch, Cer, 300x4mm, 40°dwn



- Larger 16 mm opening
- Ceramic Black Finish reduces glare
- Down 40° Angled Foot Plate
- 300 mm Shaft Length
- 4 mm bite
- Regular Footplate for Lumbar applications
- Large handle for large hands

Lrg-Hndl Punch, Cer, 300x4mm, 40°
down-angle

http://buxtonbio.com/lrgwnd_pnchcer300x4mm40dwn_1312