

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-4023

Frs-Sm Pnch, Cer, 180x2mm, 90°dwn



Photo
Pending

- Ceramic Black Finish
- 0
- Down 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Ferris-Smith Handle
- 180 mm Shaft Length

Frs-Sm Punch, Cer, 180x2mm, 90°
down-angle

http://buxtonbio.com/frssm_pnchcer180x2mm90dwn_1326