

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

847-8314

Pnch,CerR-HdlBay,180x1,40°TFP



- Ceramic
- Original
- Offset Right Handed
- Shaft Length 180mm
- Chip Set None
- 40° Angled Thin Foot Plate

Punch,Cer, Rt-Handled Bayonet,
180x1mm,40°TFP

http://buxtonbio.com/pnchcerrhdlbay180x140tfp_1733