

# Buxton BioMedical



Buxton BioMedical  
11 Melanie Lane  
East Hanover, NJ 07936  
USA

Phone: 973-560-4848  
Fax: 973-560-4999

## 848-3021

Punch, Cer, 180x2mm, 90°up, Ejec



Photo  
Pending

- Ceramic Black Finish
- 0
- Up 90° Angled Foot Plate
- Foot Plate Thickness Regular
- Chip Handling Ejector
- Small Handle
- 180 mm Shaft Length

Sm-Hndl Punch, Cer, 180x2mm, 90°  
up-ang, Ejec

[http://buxtonbio.com/punchcer180x2mm90upejec\\_1272](http://buxtonbio.com/punchcer180x2mm90upejec_1272)