

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-5261

Pnch, Cer, 200x6, 90°up, Ejec, TFP



- Ceramic Black Finish
- 0
- Up 90° Angled Foot Plate
- Foot Plate Thickness Thin
- Chip Handling Ejector
- 200 mm Shaft Length

Sm-Hndl Punch, Cer, 200x6mm, 90°up,
Ejec, TFP

http://buxtonbio.com/pnchcer200x690upejectfp_1292