

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-7631

Bay-Ker, Up, RR, Cer90°, 150x3mm



Photo
Pending

- Ceramic
- Original
- Offset Up
- Shaft Length 150mm
- Chip Set Ejector
- 40° Angled Thin Foot Plate

Bay-Ker, Up, RR, Cer90°, 150x3mm

http://buxtonbio.com/baykeruprrcer90150x3mm_1672