

Buxton BioMedical



Buxton BioMedical
11 Melanie Lane
East Hanover, NJ 07936
USA

Phone: 973-560-4848
Fax: 973-560-4999

848-7664

Bay-Ker, Up, RR, Cer40°, 150x6mm, tfp



- Ceramic
- Original
- Offset Up
- Shaft Length 150mm
- Chip Set Ejector
- 90° Angled Foot Plate

Bay-Ker, Up, RR, Cer40°, 150x6mm, tfp

http://buxtonbio.com/baykeruprrcer40150x6mmtfp_1674